SKINVIVE™
BY JUVÉDERM®

2 x 1mL

Allergan
Aesthetics

20089473
Revision 2024-11-25

ALLERGAN
Route de Promery
Zone Artisanale de Pré-Mairy
Pringy 74370 Annecy
France
Tel: +33 (0) 4 50 27 27 03

0344

20089473

en Only for professional use

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




COMPOSITION

Hyaluronic Acid gel	12.0 mg
Lidocaine hydrochloride	3.0 mg
Phosphate buffer pH 7.2 q.s.	1 mL (or g)

One syringe contains 1 mL of **SKINVIVE™ by Juvéderm®**.

DESCRIPTION

SKINVIVE™ by Juvéderm® is a sterile pyrogen-free physiological solution of cross-linked hyaluronic acid which is not of animal origin. The gel is presented in a graduated, pre-filled, disposable syringe. Each box contains two 1 mL **SKINVIVE™ by Juvéderm®** syringes, 4 single-use 32G1/2" sterile needles to be used only for injecting **SKINVIVE™ by Juvéderm®**, an instruction leaflet and a set of labels in order to ensure traceability.

STERILISATION

The contents of the **SKINVIVE™ by Juvéderm®** syringes is sterilised by moist heat. The 32G1/2" CE-marked needles are sterilised by radiation.

INDICATIONS

- **SKINVIVE™ by Juvéderm®** is an injectable implant used for the treatment, by filling, of superficial cutaneous depressions such as fine lines and for additional improvement of skin quality attributes such as hydration and elasticity.
- **SKINVIVE™ by Juvéderm®** is intended for intra-dermal injection by an authorized medical practitioner in accordance with local applicable regulation.
- The presence of lidocaine is meant to reduce the patient's pain during treatment.

CONTRA-INDICATIONS

- Do not inject **SKINVIVE™ by Juvéderm®** in the eyelids.
- Do not inject into the blood vessels (intravascular). Intravascular injection may lead to embolization, occlusion of the vessels, ischemia or infarction.
- Do not overcorrect.
- **SKINVIVE™ by Juvéderm®** must not be used in:
 - Patients suffering from untreated epilepsy;
 - Patients who tend to develop hypertrophic scarring;
 - Patients with known hypersensitivity to hyaluronic acid and/or to gram positive bacterial proteins as hyaluronic acid is produced by *Streptococcus* type bacteria;
 - Patients with known hypersensitivity to lidocaine or to amide-type local anaesthetics;
 - Patients suffering from porphyria;
 - Women who are pregnant or breastfeeding;
 - Children.
- **SKINVIVE™ by Juvéderm®** must not be used in areas presenting cutaneous inflammatory and/or infectious processes (acne, herpes, etc.).
- **SKINVIVE™ by Juvéderm®** should not be used simultaneously with laser treatment, deep chemical peels or dermabrasion. For surface peels, it is recommended not to inject the product if the inflammatory reaction generated is significant.

PRECAUTIONS FOR USE

- **SKINVIVE™ by Juvéderm®** is indicated only for intra-dermal injections.
- Medical practitioners must take into account the fact that this product contains lidocaine.
- **SKINVIVE™ by Juvéderm®** is not intended for use in breast augmentation/reconstruction.
- As a matter of general principle, injection of a medical device is associated with a risk of infection. Standard precautions associated with injectable materials shall be followed.
- There is no available clinical data about injection of **SKINVIVE™ by Juvéderm®** into an area which has already been treated with a non-ALLERGAN dermal filler.

It is recommended not to inject into a site which has been treated with a permanent implant.

No clinical data is available regarding the efficiency and tolerance of **SKINVIVE™ by Juvéderm®** injections in patients having a history of, or currently suffering from, autoimmune disease or autoimmune deficiency or being under immunosuppressive therapy. The medical practitioner shall therefore decide on the indication on a case-by-case basis, according to the nature of the disease and its corresponding treatment, and shall also ensure the specific monitoring of these patients. In particular, it is recommended that these patients undergo a preliminary skin testing for hypersensitivity, and to refrain from injecting the product if the disease is active.

There is no available clinical data concerning the tolerance of **SKINVIVE™ by Juvéderm®** injection in patients presenting a history of severe and/or multiple allergies. The medical practitioner shall therefore decide on the indication on a case-by-case basis, according to the nature of the allergy, and shall also ensure the specific monitoring of these at-risk patients. In particular, the decision may be taken to propose a skin testing for hypersensitivity or suitable preventive treatment prior to any injection. In case of history of anaphylactic shock, it is recommended not to inject the product.

Patients showing a history of streptococcal disease (recurrent sore throats, acute rheumatic fever) shall be subjected to a skin testing for hypersensitivity before any injection is administered. In the event of acute rheumatic fever with heart complications, it is recommended not to inject the product.

Patients on anti-coagulation medication or using substances that can prolong bleeding (warfarin, acetylsalicylic acid, nonsteroidal anti-inflammatory drugs, or other substances known to increase coagulation time such as herbal supplements with garlic or ginkgo biloba, etc.) must be warned of the potential increased risks of bleeding and haematomas during injection.

There is no data available regarding the safety of injecting greater amount than 20 mL of ALLERGAN dermal fillers per 60 kg (130 lbs) body mass per year.

Due to presence of lidocaine, the combination of **SKINVIVE™ by Juvéderm®** with certain drugs that reduce or inhibit hepatic metabolism (cimetidine, beta-blockers, etc.) is not recommended.

Due to presence of lidocaine, **SKINVIVE™ by Juvéderm®** should be used with caution in patients showing symptoms of cardiac conduction disorders.

Please recommend that the patient not use any makeup during the 12 hours following the injection treatment and that any extended exposure to the sun, UV rays and temperatures below 0°C be avoided, as well as any sauna or hammam sessions during the two weeks following the injection treatment.

The composition of this product is compatible with fields used for magnetic resonance imaging.

INCOMPATIBILITIES

Hyaluronic acid is known to be incompatible with quaternary ammonium salts such as benzalkonium chloride. **SKINVIVE™ by Juvéderm®** should therefore never be placed in contact with these substances or with medical-surgical instrumentation which has been treated with this type of substance.

There is no known interaction with other local anaesthetics.

UNDESIRABLE EFFECTS

The patients must be informed that there are potential side effects associated with implantation of this product, which may occur immediately or may be delayed. These include, but are not limited to:

Inflammatory reactions (redness, oedema, erythema, etc.) which may be associated with itching and/or pain on pressure and/or paresthesia, occurring after the injection. These reactions may last for a week.

Haematomas.

Induration or nodules at the injection site.

Staining or discolouration of the injection site might be observed, especially when HA dermal filler is injected too superficially and/or in thin skin (Tyndall effect).

- If immediate blanching occurs at any time during the injection, the injection should be stopped and appropriate action taken such as massaging the area until its return to a normal color.
- The degree and duration of the correction depend on the character of the defect treated, the tissue stress at the implant site, the depth of the implant in the tissue and the injection technique. The amount injected will depend on the areas which are to be corrected based on the experience of the medical practitioner.
- Do not overcorrect as injection of an excessive volume can be at the origin of some side effects such as tissue necrosis and oedema.
- A touch up (for achieving optimal correction) and/or a repeat (for maintaining optimal correction) treatment with **SKINVIVE™ by Juvéderm®** might be required.
- It is recommended to wait until side effects are resolved (with a minimal interval of 2 weeks) between two injections.
- It is important to massage the area treated after the injection in order to ensure that the substance has been uniformly distributed.

WARNINGS

- Check the expiry date on the product label.
- In the event that the content of a syringe shows signs of separation and/or appears cloudy, do not use the syringe.
- Do not re-use. Sterility of this device cannot be guaranteed if the device is re-used.
- Do not re-sterilise.
- For the needles :
Used needles must be thrown away in the appropriate containers. Do the same for the syringes. Please consult the current applicable directives to ensure their correct elimination.
- Never try to straighten a bent needle; throw it away and replace it.

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