VOLBELLA®
WITH LIDOCAINE

2 x 1mL

 **Allergan.**

73363JS10
Revision 2020-07-10

For the syringe: 

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



COMPOSITION

Hyaluronic Acid gel	15 mg
Lidocaine hydrochloride	3 mg
Phosphate buffer pH 7.2 q.s.	1 mL

One syringe contains 1 mL of **Juvéderm® VOLBELLA® with Lidocaine**.

DESCRIPTION

Juvéderm® VOLBELLA® with Lidocaine is a sterile pyrogen-free physiological solution of cross-linked hyaluronic acid which is not of animal origin. The gel is presented in a graduated, pre-filled, disposable syringe. Each box contains two 1 mL **Juvéderm® VOLBELLA® with Lidocaine** syringes, 4 single-use 32G1/2" sterile needles to be used only for injecting **Juvéderm® VOLBELLA® with Lidocaine**, an instruction leaflet and a set of labels in order to ensure traceability.

STERILISATION

The contents of the **Juvéderm® VOLBELLA® with Lidocaine** syringes is sterilised by moist heat.

The 32G1/2" needles are sterilised by radiation.

INDICATIONS

- **Juvéderm® VOLBELLA® with Lidocaine** is an injectable implant used for the treatment of any fine lines and medium-sized skin depressions due to conditions such as premature aging.
- **Juvéderm® VOLBELLA® with Lidocaine** can also be used for enhancement and pouting of the lips to correct structural defects such as asymmetry, contour deformities, volume loss...
- **Juvéderm® VOLBELLA® with Lidocaine** is intended to be used via superficial or mid-dermis injection or lips mucosa injection by an authorized medical practitioner. For the treatment of the infraorbital skin depressions (also known as tear trough area), it is recommended to inject in the submuscular / pre- periosteal plane.
- The presence of lidocaine is meant to reduce the patient's pain during treatment.

CONTRA-INDICATIONS

- Do not inject **Juvéderm® VOLBELLA® with Lidocaine** into the eyelids.
- Do not inject into the blood vessels (intravascular). Intravascular injection may lead to embolization, occlusion of the vessels, ischemia or infarction.
- Do not overcorrect.
- **Juvéderm® VOLBELLA® with Lidocaine** must not be used in:
 - Patients suffering from untreated epilepsy;
 - Patients who tend to develop hypertrophic scarring;
 - Patients with known hypersensitivity to hyaluronic acid and/or to gram positive bacterial proteins as hyaluronic acid is produced by *Streptococcus* type bacteria;
 - Patients with known hypersensitivity to lidocaine or to amide-type local anaesthetics;
 - Patients suffering from porphyria;
 - Women who are pregnant or breastfeeding;
 - Children (under 18 years).
- **Juvéderm® VOLBELLA® with Lidocaine** must not be used in areas presenting cutaneous inflammatory and/or infectious processes (acne, herpes, etc.).
- **Juvéderm® VOLBELLA® with Lidocaine** should not be used simultaneously with laser treatment, deep chemical peels or dermabrasion. For surface peels, it is recommended not to inject the product if the inflammatory reaction generated is significant.

PRECAUTIONS FOR USE

- **Juvéderm® VOLBELLA® with Lidocaine** is indicated for intradermal injection for fine lines and moderate skin depressions, injection into the mucous membrane for the lips, and submuscular / pre-periosteal injection for the tear trough/infraorbital area.
- Medical practitioners must take into account the fact that this product

contains lidocaine.

- **Juvéderm® VOLBELLA® with Lidocaine** is not intended for use in breast augmentation/reconstruction.
- As a matter of general principle, injection of a medical device is associated with a risk of infection. Standard precautions associated with injectable materials shall be followed.
- There is no available clinical data about injection of **Juvéderm® VOLBELLA® with Lidocaine** into an area which has already been treated with a non-ALLERGAN dermal filler.
- It is recommended not to inject into a site which has been treated with a permanent implant.
- No clinical data is available regarding the efficiency and tolerance of **Juvéderm® VOLBELLA® with Lidocaine** injections in patients having a history of, or currently suffering from, autoimmune disease or autoimmune deficiency or being under immunosuppressive therapy. The medical practitioner shall therefore decide on the indication on a case-by-case basis, according to the nature of the disease and its corresponding treatment, and shall also ensure the specific monitoring of these patients. In particular, it is recommended that these patients undergo a preliminary skin testing for hypersensitivity, and to refrain from injecting the product if the disease is active.
- There is no available clinical data concerning the tolerance of **Juvéderm® VOLBELLA® with Lidocaine** injection in patients presenting a history of severe and/or multiple allergies. The medical practitioner shall therefore decide on the indication on a case-by-case basis, according to the nature of the allergy, and shall also ensure the specific monitoring of these at-risk patients. In particular, the decision may be taken to propose a skin testing for hypersensitivity or suitable preventive treatment prior to any injection. In case of history of anaphylactic shock, it is recommended not to inject the product.
- No clinical data is available concerning treatment in the infraorbital region with **Juvéderm® VOLBELLA® with Lidocaine** in patients with a pre-existing tendency for infraorbital oedema. The medical practitioner should therefore consider the patient's medical history and infraorbital anatomy and physiology in evaluating whether the patient is a candidate for treatment. These patients must be warned of the potential increased risk of prolonged oedema that may occur following infraorbital treatment.
- Patients showing a history of streptococcal disease (recurrent sore throats, acute rheumatic fever) shall be subjected to a skin testing for hypersensitivity before any injection is administered. In the event of acute rheumatic fever with heart complications, it is recommended not to inject the product.
- Patients on anti-coagulation medication or using substances that can prolong bleeding (warfarin, acetylsalicylic acid, nonsteroidal anti-inflammatory drugs, or other substances known to increase coagulation time such as herbal supplements with garlic or ginkgo biloba, etc.) must be warned of the potential increased risks of bleeding and haematomas during injection.
- There is no data available regarding the safety of injecting greater amount than 20 mL of ALLERGAN dermal fillers per 60 kg (130 lbs) body mass per year.
- Due to presence of lidocaine, the combination of **Juvéderm® VOLBELLA® with Lidocaine** with certain drugs that reduce or inhibit hepatic metabolism (cimetidine, beta-blockers, etc.) is not recommended.
- Due to presence of lidocaine, **Juvéderm® VOLBELLA® with Lidocaine** should be used with caution in patients showing symptoms of cardiac conduction disorders.
- Please recommend that the patient not use any makeup during the 12 hours following the injection treatment and that any extended exposure to the sun, UV rays and temperatures below 0°C be avoided, as well as any sauna or hammam sessions during the two weeks following the injection treatment.
- The composition of this product is compatible with fields used for magnetic resonance imaging.

INCOMPATIBILITIES

Hyaluronic acid is known to be incompatible with quaternary ammonium salts such as benzalkonium chloride. **Juvéderm® VOLBELLA® with Lidocaine** should therefore never be placed in contact with these substances or with

medical-surgical instrumentation which has been treated with this type of substance.

There is no known interaction with other local anaesthetics.

UNDESIRABLE EFFECTS

The patients must be informed that there are potential side effects associated with implantation of this product, which may occur immediately or may be delayed. These include, but are not limited to:

- Inflammatory reactions (redness, oedema, erythema, etc.) which may be associated with itching, pain on pressure, and/or paresthesia, occurring after the injection. These reactions may last for a week. In particular, injection in the mucous membrane of the lips or infraorbital area may cause more oedema and bruising due to the specific physiology of these tissues. A preventive anti-inflammatory treatment by a medical practitioner can be recommended.
- Haematomas.
- Induration or nodules at the injection site.
- Staining or discolouration of the injection site might be observed, especially when HA dermal filler is injected too superficially and/or in thin skin (Tyndall effect).
- Poor effect or weak filling effect.
- Rare but serious adverse events associated with intravascular injection of dermal fillers in the face and tissue compression have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis and damage to underlying structures. Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in the vision, signs of stroke, blanching of the skin or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate medical practitioner specialist should an intravascular injection occur. Abscesses, granuloma and immediate or delayed hypersensitivity after hyaluronic acid and/or lidocaine injections have also been reported. It is therefore advisable to take these potential risks into account.
- Patients must report inflammatory reactions which persist for more than one week, or any other side effect which develops, to their medical practitioner as soon as possible. The medical practitioner should use an appropriate treatment.
- Any other undesirable side effects associated with injection of **Juvéderm® VOLBELLA® with Lidocaine** must be reported to the distributor and/or to the manufacturer.

METHOD OF USE – POSOLOGY

- This product is designed to be injected into the dermis, in the mucous membrane of the lips, or in the submuscular / pre-periosteal plane in the infraorbital area, by an authorized medical practitioner in accordance with local, applicable regulation(s). In order to minimize the risks of potential complications, and as precision is essential to a successful treatment, the product should only be used by medical practitioners who have experience in injection techniques for the area of use, and by those who are knowledgeable about the anatomy and physiology at and around the site of injection.
- Use of the supplied 32G1/2" needle is recommended. However, depending on the medical practitioner's preferred injection technique, it is possible to use:
 - a 30G sterile cannula (please refer to the list hereunder). Choice of cannula length is determined by the user according to his/her injection technique,
 - a 30G1/2" sterile needle

Material Number	Description
94323/ HPC30019ACSH	Easyflow System-20* cannula 30G x 19mm.
94324/ HPC30025ACSH	Easyflow System-20* cannula 30G x 25mm.
HPC-30013A	30G1/2" TSK Hypodermic needle

- Contra-Indications, Method of use, Precautions for use and Warnings defined for the needle in this leaflet apply also to the cannula referenced

above if used with this product.

• **Juvéderm® VOLBELLA® with Lidocaine** is to be used as supplied. Modification or use of the product outside the Directions for Use may adversely impact the sterility, homogeneity and performance of the product and it can therefore no longer be assured.

• Prior to treatment, medical practitioners shall inform their patients about the product's indications, contra-indications, incompatibilities and potential undesirable effects/risks associated with dermal fillers injection and ensure that patients are aware of signs and symptoms of potential complications.

• The area to be treated should be disinfected thoroughly prior to the injection.

• Remove tip cap by pulling it straight off the syringe as shown in fig. 1. Then firmly push the needle provided in the box (fig. 2) into the syringe, screwing it gently clockwise. Twist once more until it is fully locked and has the needle cap in the position shown in fig. 3. If the needle cap is positioned as shown in fig. 4, it is incorrectly attached. Next, remove the protective cap by holding the body of the syringe in one hand, the protective cap in the other, as shown in fig. 5, and pulling the two hands in opposite directions. Prior to injecting, depress the plunger rod until the product flows out of the needle.

Inject slowly and apply the least amount of pressure necessary.

If the needle is blocked, do not increase the pressure on the plunger rod. Instead, stop the injection and replace the needle.

Failure to comply with these precautions could cause a disengagement of the needle and/or product leakage at luer-lock level and/or increase the risk of vascular compromise.

• After needle insertion and before injection, it is recommended to withdraw slightly the plunger to aspirate and verify the needle is not intravascular.

• For injections in the infraorbital area, the following technique is recommended:

- Insert the needle below the orbital rim perpendicular to the skin surface and advance until contacting the periosteum

- Slowly inject in the submuscular/pre-periosteal plane along the infraorbital hollow with the needle using a depot injection technique. Reposition the needle as needed in between depot injections, always ensuring that the needle is directly below the muscle and aspirating prior to injection to ensure the needle tip is not within a vessel.

- If a cannula is used instead of a needle, use the introducer needle to create an insertion point for the cannula below the orbital rim and inject using a retrograde tunneling technique.

• If immediate blanching occurs at any time during the injection, the injection should be stopped and appropriate action taken such as massaging the area until its return to a normal color.

• The degree and duration of the correction depend on the character of the defect treated, the tissue stress at the implant site, the depth of the implant in the tissue and the injection technique. The amount injected will depend on the areas which are to be corrected based on the experience of the medical practitioner.

• Do not overcorrect as injection of an excessive volume can be at the origin of some side effects such as tissue necrosis and oedema.

• A touch up (for achieving optimal correction) and/or a repeat (for maintaining optimal correction) treatment with **Juvéderm® VOLBELLA® with Lidocaine** might be required.

• It is recommended to wait until side effects are resolved (with a minimal interval of 2 weeks) between two injections.

• It is important to massage the area treated after the injection in order to ensure that the substance has been uniformly distributed.

WARNINGS

• Check the expiry date on the product label.

• In the event that the content of a syringe shows signs of separation and/or appears cloudy, do not use the syringe.

• Do not re-use. Sterility of this device cannot be guaranteed if the device is re-used.

• Do not re-sterilise.

• For the needles :

- Used needles must be thrown away in the appropriate containers. Do the

same for the syringes. Please consult the current applicable directives to ensure their correct elimination.

- Never try to straighten a bent needle; throw it away and replace it.

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